

# Top-down Proteomics: Comprehensive Analysis



10-80 kDa

## Intact Protein

Complete sequence + PTMs

VS



Bottom-up

## Digested Peptides

Fragmented before analysis

### Intact Protein Analysis

- No digestion required
- Analyze whole proteins
- 10-80 kDa typical range

### PTM Preservation

- Complete modification pattern
- Combinatorial PTM analysis
- Proteoform characterization

### Technical Challenges

- Requires high resolution
- Complex spectra interpretation

### Native MS


- Preserve non-covalent interactions
- Protein complexes

- Lower sensitivity than bottom-up

- Quaternary structure information

## Detailed Analysis

### 1 Intact Protein Analysis



#### Workflow Advantages

Top-down proteomics analyzes intact proteins without enzymatic digestion, preserving the complete molecular context. This approach maintains connectivity between amino acids and modifications throughout the entire protein sequence.

#### Mass Range Considerations

The effective mass range for top-down proteomics is typically 10-80 kDa, with optimal performance between 20-50 kDa. Larger proteins present challenges in ionization efficiency and spectral complexity.

**Sample Preparation:** Minimal preparation reduces analysis time and sample loss

**Direct Analysis:** Proteins are ionized and analyzed in their intact form

**Complete Information:** Full sequence coverage with positional PTM information

**Small Proteins (10-20 kDa):** Excellent sensitivity and resolution

**Medium Proteins (20-50 kDa):** Optimal balance of coverage and signal

**Large Proteins (50-80 kDa):** Challenging but achievable with advanced instrumentation

### Instrumentation Requirements

High-resolution mass spectrometers are essential for resolving the complex charge state distributions and isotope patterns of intact proteins.

**FT-ICR MS:** Ultra-high resolution (>500,000) for large proteins

**Orbitrap:** High resolution (240,000+) with good sensitivity

**Q-TOF:** Lower resolution but faster acquisition

### Key Applications

Intact protein analysis is particularly valuable for characterizing protein variants, isoforms, and complex modification patterns that would be lost in bottom-up approaches.

**Antibody Characterization:** Complete mAb analysis including glycosylation

**Biomarker Discovery:** Disease-specific protein variants



**Quality Control:** Biopharmaceutical product consistency

⚡ **Key Advantage:** Eliminates inference problems associated with peptide-level analysis

2

## PTM Preservation & Proteoform Analysis

## Proteoform Identification




### Complete PTM Patterns

Top-down proteomics uniquely preserves the complete PTM landscape of each protein molecule, revealing combinations of modifications that co-exist on individual proteoforms.

**Combinatorial PTMs:** Identifies which modifications occur together

**Stoichiometry:** Determines relative abundance of each proteoform

**Positional Information:** Precise localization of modifications

 **Bottom-up approaches cannot determine which PTMs occur on the same molecule**

### Proteoform Characterization

A proteoform is a specific molecular form of a protein arising from genetic variations, alternative splicing, or PTMs. Top-down MS is the gold standard for proteoform identification.

**Isoforms:** Alternative splice variants and sequence variants

**Modified Forms:** Unique PTM combinations

**Processed Forms:** Cleavage products and mature proteins

### Common PTMs Detected

### Biological Significance

Top-down proteomics excels at detecting and localizing diverse post-translational modifications across the entire protein sequence.

**Phosphorylation (+80 Da):** Signaling pathway regulation

**Acetylation (+42 Da):** Transcriptional regulation, histone marks

**Methylation (+14 Da):** Chromatin regulation, protein stability

**Ubiquitination (+8.5 kDa):** Protein degradation signals

**Glycosylation (variable):** Protein folding, stability, function


Understanding proteoform diversity is crucial because different proteoforms can have distinct biological functions, localizations, or disease associations.

**Cell Signaling:** Phosphorylation cascades and crosstalk

**Epigenetics:** Histone code readout



**Disease Markers:** Aberrant PTM patterns in disease

**Drug Targets:** Therapeutic intervention points

 **A single gene can produce 10-100+ distinct proteoforms through PTMs**

### 3 Technical Challenges & Solutions

#### Resolution & Sensitivity Requirements




## Resolution Requirements

Intact proteins produce complex isotope patterns and closely spaced charge states that demand ultra-high resolution to resolve accurately.

**Isotope Resolution:**  $R > 100,000$  needed for 50 kDa proteins

**Charge State Separation:** Distinguish neighboring charge states

**Mass Accuracy:**  $<5$  ppm for confident identification

 **Insufficient resolution leads to ambiguous proteoform assignments**

## Sensitivity Limitations

Top-down MS typically requires 100-1000× more sample than bottom-up approaches due to lower ionization efficiency and signal dispersion across charge states.

**Ionization Efficiency:** Decreases with increasing protein size

**Signal Dilution:** Ion current distributed across many charge states

**Dynamic Range:** Challenges detecting low-abundance proteoforms

## Spectral Deconvolution

Large proteins exhibit multiple charge states (typically 20-50+), creating overlapping spectral features that require sophisticated algorithms to interpret.

**Charge State Distribution:** Gaussian-like patterns centered at high  $z$

**Data Processing:** Automated deconvolution to zero-charge mass

**Software Tools:** ProSight, TDPortal, Informed-Proteomics

## Overcoming Challenges


Recent technological advances have significantly improved the feasibility and throughput of top-down proteomics.

**Instrument Advances:** Higher sensitivity Orbitraps, 21T FT-ICR



**Front-end Separation:** RPLC, CZE, GELFrEE

**Enrichment Strategies:** Immunoprecipitation, affinity capture

**Novel Fragmentation:** UVPD, AI-ETD for improved coverage

 **Combining multiple separation dimensions improves coverage by 10-fold**

## Preserving Native Structure



## Native MS Principles

Native mass spectrometry maintains proteins in near-physiological conditions during ionization, preserving non-covalent interactions and native conformations in the gas phase.

**Buffer Conditions:** Aqueous solutions with physiological pH and salt

**Gentle Ionization:** Nano-ESI with minimal desolvation energy

**Structure Preservation:** Complexes remain intact during transfer

## Protein Complex Analysis

Native MS excels at characterizing multi-protein complexes, providing information on stoichiometry, assembly state, and stability that is difficult to obtain by other methods.

**Subunit Composition:** Identify all complex components

**Stoichiometry:** Determine exact subunit ratios

**Assembly States:** Detect oligomers and aggregates

**Heterogeneity:** Reveal co-existing assembly forms



Native MS bridges structural biology and mass spectrometry

## Quaternary Structure Information

Native MS provides unique insights into the architecture and dynamics of protein complexes at the quaternary structure level.

**Complex Size:** MDa-scale assemblies (ribosomes, proteasomes)

**Binding Partners:** Protein-protein, protein-ligand interactions

**Conformational States:** Distinguish active vs. inactive forms

**Stability Assessment:** Measure dissociation constants (Kd)

## Applications & Examples

Native MS has become indispensable in structural biology and biopharmaceutical development for understanding protein interactions and assembly.

**Antibody-Antigen:** Characterize mAb binding and valency

**Viral Capsids:** Assembly mechanisms and stability

**Membrane Proteins:** Lipid and detergent interactions

**Drug Discovery:** Screen ligand binding and selectivity

**Quality Control:** AAV, VLP, and vaccine characterization



Native MS is FDA-recognized for biotherapeutic characterization