



# Mass Spectrometry Basics



## Ion Generation

ESI / MALDI



## Mass Analysis

$m/z$  Separation



## Detection

Signal Processing

### Resolution & Accuracy



- High resolution: distinguish similar masses
- Mass accuracy: parts per million (ppm)
- Critical for peptide identification

### Scan Modes



- Full scan: entire mass range
- Selected ion monitoring (SIM)
- Data-dependent acquisition (DDA)

### Sensitivity



- Femtomole to attomole detection
- Dynamic range: 3-5 orders of magnitude
- Low abundance protein detection

### Applications



- Protein identification
- Quantification
- PTM analysis

# 1. Resolution & Accuracy

## High Resolution

m/z 1000.50 m/z 1000.52



$\Delta m = 0.02$  Da

## Low Resolution

Unresolved



### Resolution Formula:

$$R = m / \Delta m$$

where m = mass and  $\Delta m$  = peak width at half maximum (FWHM)

## Mass Resolution

The ability to distinguish between ions of similar mass-to-charge ratios. Higher resolution allows separation of peaks that differ by small mass units.

### Example:

Orbitrap mass analyzers achieve resolutions >100,000, enabling discrimination between peptides differing by **0.001 Da**.

## Mass Accuracy

The closeness of measured mass to the true mass, typically expressed in parts per million (ppm). Essential for confident peptide identification.

$$\text{Mass Accuracy (ppm)} = [(m_{\text{measured}} - m_{\text{theoretical}}) / m_{\text{theoretical}}] \times 10^6$$


### Typical Values:

Modern instruments: **<5 ppm** accuracy


High-resolution MS: **<1 ppm** accuracy

## 2. Scan Modes

### Full Scan (MS1)



### Selected Ion Monitoring (SIM)



### Data-Dependent Acquisition (DDA)



### Full Scan Mode

Measures all ions across a specified  $m/z$  range. Provides comprehensive overview but lower sensitivity for individual ions.

### Selected Ion Monitoring (SIM)

Focuses on specific  $m/z$  values of interest. Increases sensitivity and detection limit by 10-100 fold for target analytes.

#### Application:

Targeted quantification of known peptides in complex mixtures, such as biomarker validation studies.

### Data-Dependent Acquisition (DDA)

Automatically selects the most abundant ions from MS1 scan for fragmentation (MS2). Enables protein identification through peptide sequencing.


#### Workflow:

1. Survey scan (MS1)
2. Select top N peaks

3. Fragment selected ions (MS2)

4. Return to survey scan

### 3. Sensitivity



#### Detection Sensitivity

Modern MS instruments can detect proteins at femtomole ( $10^{-15}$  mol) to attomole ( $10^{-18}$  mol) levels, enabling analysis of rare proteins and post-translational modifications.

##### Practical Context:

**1 femtomole** = approximately 600,000 molecules

**1 attomole** = approximately 600 molecules

Sufficient for single-cell proteomics

#### Dynamic Range

The ratio between the most and least abundant detectable proteins. Typical range:  **$10^3$  to  $10^5$**  (3-5 orders of magnitude).

#### Challenges

- High-abundance proteins (e.g., albumin) can mask low-abundance proteins
- Sample depletion strategies often required
- Ion suppression effects in complex mixtures

##### Solution:

Fractionation techniques (e.g., HPLC) combined with enrichment methods improve detection of low-abundance proteins by **100-1000 fold**.

## 4. Applications in Proteomics

### 1. Protein Identification

Protein → Digest → Peptides



### 2. Quantitative Proteomics

#### Label-Free

- Spectral counting
- Peak intensity
- No labeling required

#### Isotope Labeling

- SILAC, TMT, iTRAQ
- Multiplexing
- High precision

### 3. Post-Translational Modification (PTM)

Phospho  
+80 Da

Acetyl  
+42 Da

Ubiquitin  
+114 Da

Mass shift detection enables PTM identification

### Protein Identification

Proteins are digested into peptides, which are analyzed by MS/MS. Fragment ion patterns are matched against databases for identification.

#### Database Search:

Mascot, SEQUEST, MaxQuant algorithms match experimental spectra to theoretical peptide fragments. Typical identification: **>2 unique peptides** per protein.

### Quantitative Analysis

Compares protein abundance across samples using label-free or isotope labeling approaches. Essential for biomarker discovery and pathway analysis.

### PTM Characterization

Detects modifications through characteristic mass shifts. Phosphorylation, acetylation, methylation, ubiquitination are routinely analyzed.

#### Clinical Application:

PTM analysis in cancer research: Aberrant phosphorylation patterns identify activated signaling pathways and therapeutic targets.