


# Flow Cytometry Principles

## Flow Cytometer System



Analysis Speed

**10,000/s**

Typical Lasers

**405, 488, 561, 640 nm**

Parameters/Cell

**20-50+**

## Fluidics System

Hydrodynamic focusing creates single-cell stream  
Sheath fluid (PBS) surrounds sample  
Laminar flow for precise alignment

## Laser Excitation

Multiple lasers for multicolor detection  
Common: 405, 488, 561, 640 nm  
Each excites different fluorophores

## Detection Channels

**FSC:** Forward scatter (cell size)  
**SSC:** Side scatter (granularity)  
**FL1-FLn:** Fluorescence PMTs

## Compensation

Corrects spectral overlap between fluorophores  
Single-color controls essential  
Software or hardware compensation

## Applications:

Immunophenotyping • Cell cycle analysis • Apoptosis detection • Rare cell identification • Biomarker expression