

Library Preparation Methods for RNA-Seq

PolyA Selection

Enriches mRNA by capturing poly-adenylated transcripts

Ribosomal Depletion

Removes rRNA to capture all RNA types including non-coding

Strand Specificity

Preserves information about which DNA strand was transcribed

UMI Incorporation

Unique Molecular Identifiers enable accurate quantification


3' Tag-seq

Sequences only 3' ends - cost effective for counting

Full-length Coverage

Complete transcript coverage for isoform analysis

Trade-off: **Cost vs. Information Content**



Detailed Method Explanations

1 PolyA Selection

Overview

PolyA selection is the most commonly used method for mRNA enrichment in RNA-seq. It takes advantage of the polyadenylated (polyA) tail present on most mature eukaryotic mRNAs. Oligo(dT) beads are used to bind and capture these polyA-tailed transcripts, effectively enriching for coding RNA while removing ribosomal RNA and other non-polyadenylated species.

Key Features

- ▶ Captures ~1-2% of total RNA (mostly mRNA)
- ▶ Removes ~90% of ribosomal RNA
- ▶ Standard method for gene expression studies
- ▶ Compatible with degraded samples (with limitations)

✓ Advantages

- ▶ Cost-effective
- ▶ Simple protocol

✗ Limitations

- ▶ Misses non-polyA RNAs


- ▶ High enrichment efficiency
- ▶ Well-established method

- ▶ 3' bias with degradation
- ▶ Excludes bacterial transcripts
- ▶ May lose histone mRNAs

Best Use Cases

Standard gene expression profiling, differential expression analysis, transcript quantification in high-quality eukaryotic samples, and when focusing exclusively on protein-coding genes.

Total RNA



2 Ribosomal RNA Depletion

Overview

Ribosomal depletion removes rRNA (which comprises 80-90% of total RNA) using sequence-specific probes that hybridize to ribosomal sequences. Unlike polyA selection, this method retains all other RNA types including non-polyadenylated transcripts, making it ideal for comprehensive transcriptome

analysis including bacterial samples, long non-coding RNAs, and degraded samples.

Key Features

- ▶ Removes rRNA while preserving all other RNA species
- ▶ Works with prokaryotic and eukaryotic samples
- ▶ Retains non-coding RNAs (lncRNA, circRNA, etc.)
- ▶ Better for degraded samples than polyA selection

✓ Advantages

- ▶ Captures all RNA types
- ▶ No 3' bias
- ▶ Works with bacteria
- ▶ Better for FFPE samples


✗ Limitations

- ▶ More expensive
- ▶ Incomplete rRNA removal
- ▶ Complex data analysis
- ▶ Species-specific probes

Best Use Cases

Non-coding RNA studies, bacterial transcriptomics, degraded or FFPE samples, comprehensive transcriptome profiling including intronic and intergenic regions, and when polyA selection is not suitable.

Ribosomal Depletion Process



Overview

Strand-specific (stranded) RNA-seq preserves information about which DNA strand a transcript originated from. This is crucial because genes can overlap, and antisense transcripts can regulate gene expression. The method typically uses dUTP incorporation during second-strand synthesis, which is later degraded, ensuring only the original strand is sequenced.

Key Features

- ▶ Maintains strand-of-origin information
- ▶ Resolves overlapping genes on opposite strands
- ▶ Detects antisense transcription
- ▶ Improves accuracy in transcript quantification


✓ Advantages

- ▶ Distinguishes sense/antisense
- ▶ Better annotation accuracy
- ▶ Identifies regulatory RNAs

✗ Limitations

- ▶ Slightly higher cost
- ▶ More complex protocol
- ▶ Additional processing step
- ▶ Requires compatible

Strand-Specific Protocol



- ▶ Resolves overlapping genes

analysis

Best Use Cases

Discovery of antisense transcripts, analysis of overlapping gene pairs, regulatory RNA studies, transcript annotation improvement, and any study requiring precise strand information for complex genomic regions.

4 UMI (Unique Molecular Identifier) Incorporation

Overview

UMIs are short random nucleotide sequences (typically 8-12 bp) added to each RNA molecule before amplification. During sequencing, PCR duplicates can be identified and removed by their identical UMI sequences, allowing accurate counting of original molecules. This dramatically improves quantification accuracy by distinguishing biological duplicates from PCR artifacts.

Key Features

- ▶ Random barcodes tag each individual molecule
- ▶ Enables removal of PCR duplicates

- ▶ Provides absolute molecule counting
- ▶ Essential for single-cell RNA-seq

✓ Advantages

- ▶ Accurate quantification
- ▶ Removes amplification bias
- ▶ Better for low-input samples
- ▶ Improves reproducibility

✗ Limitations

- ▶ Increases sequencing cost
- ▶ Requires specialized analysis
- ▶ May reduce mapping rate
- ▶ UMI collisions possible

Best Use Cases

Single-cell RNA-seq, low-input RNA samples, precise quantification studies, detection of rare transcripts, clinical applications requiring high accuracy, and any experiment where PCR bias is a concern.

UMI Workflow

1. Original RNA molecules



2. Add UMIs during RT



3. PCR Amplification



4. Computational deduplication



Accurate molecule count: 2 original molecules

5

3' Tag-seq (3' End Sequencing)

Overview

3' Tag-seq sequences only the 3' end of transcripts (typically 50-100 bp from the polyA tail), making it highly cost-effective for gene expression quantification. Since most genes can be identified from their 3' UTR, this method provides accurate counting at a fraction of the cost of full-length sequencing. It's ideal for differential expression studies where transcript structure information is not needed.

Key Features

- ▶ Sequences only 3' terminal fragments
- ▶ Significantly reduced sequencing costs
- ▶ Maintains accurate gene-level quantification
- ▶ Compatible with UMIs for enhanced accuracy


✓ Advantages

- ▶ Highly cost-effective
- ▶ Less sequencing depth needed
- ▶ Good for large sample sets
- ▶ Robust to 5' degradation

✗ Limitations

- ▶ No isoform information
- ▶ Cannot detect SNPs
- ▶ Limited to gene counting
- ▶ Requires good 3' annotation

3' Tag-seq Strategy



Cost Comparison

Full-length sequencing \$\$

3' Tag-seq \$

Key Benefits

- ✓ 70-80% cost reduction
- ✓ More samples per budget • ✓ Faster analysis

Best Use Cases

Large-scale differential expression studies, population-level gene expression screening, cost-constrained experiments, single-cell RNA-seq, studies focusing solely on gene abundance without requiring isoform or variant information.

6 Full-length Transcript Coverage

Overview

Full-length RNA-seq captures the entire length of transcripts from 5' to 3' end, providing comprehensive information about transcript structure, isoforms, and sequence variants. This method uses specialized protocols to minimize 3' bias and ensure uniform coverage across the entire transcript length, enabling detection of alternative splicing, fusion genes, and RNA editing events.

Key Features

- ▶ Complete transcript coverage (5' to 3')
- ▶ Detects alternative splicing and isoforms
- ▶ Identifies SNPs and RNA editing sites
- ▶ Enables allele-specific expression analysis

✓ Advantages

- ▶ Complete transcript information
- ▶ Isoform identification
- ▶ Detects fusion transcripts
- ▶ Variant calling capability

✗ Limitations

- ▶ Most expensive method
- ▶ Requires high sequencing depth
- ▶ Complex data analysis
- ▶ Needs high-quality RNA

Best Use Cases

Alternative splicing analysis, novel isoform discovery, fusion gene detection in cancer, RNA editing studies, allele-specific expression, comprehensive transcriptome annotation, and any research requiring complete transcript-level information.

Full-length Coverage Analysis

