


LOINC for Laboratory Tests

LOINC Six-Part Structure

2345-7: Glucose [Mass/volume] in Serum



Common LOINC Code Examples

2160-0

Creatinine in Serum

718-7

Hemoglobin in Blood

2951-2

Sodium in Serum

4544-3

Hematocrit in Blood



Test Categories

- Chemistry tests
- Hematology & coagulation
- Microbiology cultures
- Serology & immunology
- Molecular pathology



Panel Organization

- Basic Metabolic Panel (BMP)
- Complete Blood Count (CBC)
- Comprehensive Metabolic Panel
- Lipid panel
- Liver function tests



LOINC Properties

- MCnc: Mass concentration
- NCnc: Number concentration
- Prid: Presence/Identity



UCUM Units

- mg/dL, mmol/L (chemistry)
- 10³/uL (cell counts)
- IU/L (enzymes)

- Titr: Titer (dilution)
- Arb: Arbitrary units

- Standardized unit conversion
- Reference range mapping


Detailed Category Descriptions



Test Categories

LOINC provides comprehensive coverage of laboratory test categories, each representing different analytical domains in clinical diagnostics. These categories enable systematic organization and identification of all laboratory observations.

Laboratory Test Category Hierarchy



Real-World Example: Chemistry Test

Test: Glucose measurement in serum/plasma

LOINC Code: 2345-7

Full Name: Glucose [Mass/volume] in Serum or Plasma

Clinical Use: Diabetes screening, monitoring, diagnosis of hyperglycemia/hypoglycemia

Units: mg/dL or mmol/L (UCUM standard)

Real-World Example: Hematology Test

Test: Hemoglobin measurement in whole blood

LOINC Code: 718-7

Full Name: Hemoglobin [Mass/volume] in Blood

Clinical Use: Anemia diagnosis, blood loss assessment, monitoring chronic disease


Units: g/dL or g/L (UCUM standard)




Panel Organization

Laboratory test panels group related tests together for efficient clinical assessment. LOINC provides specific codes for both individual tests and complete panels, enabling comprehensive documentation of multi-test orders.

Common Laboratory Panels



Panel Structure in LOINC



Real-World Example: Basic Metabolic Panel (BMP)

Panel LOINC Code: 51990-0

Full Name: Basic metabolic panel - Blood

Components: 8 individual tests (Glucose, Sodium, Potassium, Chloride, CO₂, BUN, Creatinine, Calcium)

Clinical Use: Routine metabolic screening, kidney function assessment, electrolyte balance

Advantage: Single order code generates all component test orders

Real-World Example: Complete Blood Count (CBC)

Panel LOINC Code: 58410-2

Full Name: Complete blood count (hemogram) panel - Blood by Automated count

Components: WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelets, RDW

Clinical Use: Anemia diagnosis, infection screening, general health assessment

Result Format: Each component has its own LOINC code for reporting

LOINC Properties

The Property axis in LOINC defines what is being measured about the analyte. This critical component distinguishes different measurement types for the same substance, ensuring precise test identification.

LOINC Property Types

MCnc

Mass Concentration

Mass per unit volume
Most common property

Example Units:
mg/dL, g/L, mmol/L

Example Test:
Glucose: 2345-7

NCnc

Number Concentration

Count per unit volume
Used for cell counts

Example Units:
 $10^3/\mu\text{L}$, $10^6/\mu\text{L}$

Example Test:
WBC count: 6690-2

Prid

Presence/Identity

Qualitative detection
Present/absent results

Example Results:
Positive, Negative

Example Test:
Bacteria ID: 600-7

Titr

Titer (dilution)

Serial dilution testing
Antibody levels

Example Results:
1:8, 1:64, 1:256

Example Test:
ANA titer: 5048-9

Arb

Arbitrary Units

Non-standard units
Immunoassays

Example Units:
[arb'U]/mL, U

Example Test:
IgE: 19113-0

Additional Properties

- **ACnc**: Activity concentration (enzymes)
- **CCnc**: Catalytic concentration
- **MFr**: Mass fraction (percentage)
- **NFr**: Number fraction
- **Rto**: Ratio
- **Vol**: Volume
- **Time**: Measurement of time

Real-World Example: Mass Concentration (MCnc)

Test: Glucose in Serum

LOINC Code: 2345-7

Property: MCnc (Mass Concentration)

Measurement: Mass of glucose per volume of serum

Units: mg/dL (conventional) or mmol/L (SI units)

Calculation: If result is 100 mg/dL, this means 100 milligrams of glucose per deciliter of serum

Real-World Example: Number Concentration (NCnc)

Test: White Blood Cell Count

LOINC Code: 6690-2

Property: NCnc (Number Concentration)

Measurement: Number of WBCs per unit volume

Units: $10^3/\mu\text{L}$ (thousands per microliter) or $10^9/\text{L}$


Interpretation: Result of 7.5 means 7,500 white blood cells per microliter



UCUM Units (Unified Code for Units of Measure)

UCUM provides a standardized system for expressing units of measure in clinical laboratory testing. Integration with LOINC ensures consistent, unambiguous reporting of test results across different healthcare systems.

UCUM Unit Categories



Unit Conversion Example: Glucose



Real-World Example: Chemistry Test Units

Test: Serum Creatinine

LOINC Code: 2160-0

Conventional Units: mg/dL (milligrams per deciliter)

SI Units: $\mu\text{mol/L}$ (micromoles per liter)

Conversion Factor: $\text{mg/dL} \times 88.4 = \mu\text{mol/L}$

Example: $1.2 \text{ mg/dL} = 106 \mu\text{mol/L}$

Reference Range: 0.7-1.3 mg/dL (62-115 $\mu\text{mol/L}$) for males

Real-World Example: Hematology Cell Count Units

Test: White Blood Cell Count

LOINC Code: 6690-2

UCUM Units: $10^3/\mu\text{L}$ (thousands per microliter)

Alternative: $10^9/\text{L}$ (billions per liter)

Conversion: $1 \times 10^3/\mu\text{L} = 1 \times 10^9/\text{L}$

Example Result: $7.5 \times 10^3/\mu\text{L}$ means 7,500 cells per microliter

Reference Range: $4.5\text{--}11.0 \times 10^3/\mu\text{L}$

Real-World Example: Enzyme Activity Units

Test: Alanine Aminotransferase (ALT)

LOINC Code: 1742-6

UCUM Units: IU/L or U/L (International Units per liter)

SI Alternative: ukat/L (microkatal per liter)

Conversion: $1 \text{ U/L} = 0.0167 \text{ ukat/L}$

Temperature: Measured at 37°C (body temperature)

Reference Range: 7-56 U/L (varies by laboratory)

LOINC® is maintained by the Regenstrief Institute and is freely available for use.

For complete LOINC database and documentation, visit **loinc.org**