


CITE-seq

Cellular Indexing of Transcriptomes and Epitopes by Sequencing

CITE-seq Workflow Overview



Key Benefits of CITE-seq

- | | |
|--|---|
| 1 Enhanced Cell Type Resolution Protein markers provide clearer cell identity | 4 Low-abundance Transcripts Proteins detect markers missed by RNA-seq |
| 2 Validation of Gene Expression Direct protein measurement validates RNA data | 5 Platform Compatibility Works with standard 10X Genomics workflow |
| 3 Post-translational Information Captures protein modifications not visible in RNA | 6 Same Cell Barcode RNA and protein data share cell identity |

💡 CITE-seq bridges the gap between transcriptomics and proteomics at single-cell resolution

1. Antibody-Derived Tags (ADT) Technology

CITE-seq's core innovation lies in the use of oligonucleotide-conjugated antibodies, known as Antibody-Derived Tags (ADTs). These specialized reagents enable simultaneous measurement of surface protein expression alongside gene expression in individual cells.

ADT Structure and Design

ADT Molecular Architecture



- **Antibody:** Recognizes specific cell surface protein (e.g., CD3, CD14, CD19)
- **Linker:** Stable chemical bond connecting antibody to DNA (often via streptavidin-biotin)

- **ADT Barcode:** Unique sequence identifying the antibody (12-15 bp)
- **Poly-A Tail:** Enables capture by oligo-dT primers on 10X beads

How ADTs Work

ADTs function by combining the specificity of antibody-antigen recognition with the quantitative power of DNA sequencing. When cells are stained with a cocktail of ADT-conjugated antibodies, each antibody binds to its cognate surface protein. During droplet encapsulation and reverse transcription, the poly-A tail of the ADT oligonucleotide is captured by the same oligo-dT primers that capture mRNA, ensuring that protein and RNA measurements share the same cell barcode.

Designing an Effective ADT Panel

Marker Selection

- Choose well-established surface markers
- Include lineage-defining proteins
- Target stable, abundant proteins

Antibody Quality

- Validate specificity by flow cytometry
- Use clone-validated antibodies
- Avoid cross-reactive antibodies

Panel Size

- Start with 20-50 markers for immune cells
- Scale up to 100+ for comprehensive profiling
- Consider cost vs. information gain

Barcode Design

- Ensure sufficient Hamming distance
- Avoid homopolymers (AAAA, GGGG)
- Balance GC content (~50%)

Staining Optimization

- Titrate antibody concentrations
- Stain at 4°C for 30 minutes
- Wash thoroughly to remove excess

Controls

- Include isotype controls
- Use unstained samples
- Add cell hashing for multiplexing


Key Advantages of ADT Technology

- **Direct Protein Measurement:** Captures actual protein expression, not just mRNA abundance
- **Higher Sensitivity:** Antibodies can detect low-abundance proteins that may not be captured in RNA data
- **Same Cell Resolution:** Both RNA and protein data are linked by the same cell barcode
- **Reduced Dropout:** Protein measurements are less affected by technical dropout compared to sparse scRNA-seq data
- **Functional States:** Captures activation markers and post-translational modifications

2. Technical Workflow and Implementation

CITE-seq integrates seamlessly with standard droplet-based single-cell RNA-seq platforms, particularly the 10X Genomics Chromium system. The workflow requires minimal modifications to existing protocols while adding a powerful protein measurement dimension.

Step-by-Step Protocol



Critical Technical Considerations

- **ADT Library Optimization:** The ADT library is much smaller and simpler than the cDNA library, requiring careful titration to avoid over-amplification. Typically, ADT libraries are sequenced with 5-10% of total sequencing reads.

- **Size Selection:** Proper size selection is crucial to separate ADT libraries (120-180 bp) from cDNA libraries (300-600 bp). SPRI bead cleanup at 0.6x ratio effectively separates these populations.
- **Washing Stringency:** Thorough washing after antibody staining is essential to remove unbound ADTs, which can create high background signal. Three washes with PBS + 1% BSA are recommended.
- **Cell Concentration:** Maintain cell concentration at recommended levels (700-1,200 cells/ μ L) to ensure proper encapsulation efficiency and minimize doublets.

Sequencing and Data Output

Sequencing Read Structure

Gene Expression Read

Read 1 (28bp)

Cell Barcode + UMI

Read 2 (91bp)

cDNA insert (gene sequence)

Example: ACGTACGT...GCTA (maps to GAPDH)

ADT Read

Read 1 (28bp)

Cell Barcode + UMI

Read 2 (15bp)

ADT barcode

Example: ACGTGCATCGATCGA (CD3 antibody)

Sequencing Depth Recommendations

Gene Expression Library

Recommended depth: 20,000-50,000 reads per cell

ADT Library

Recommended depth: 1,000-5,000 reads per cell

Purpose: Capture transcriptome diversity
Typical output: 2,000-5,000 genes per cell


Purpose: Quantify protein markers
Typical output: 20-200 proteins per cell

3. Applications and Multi-Modal Data Integration

CITE-seq enables sophisticated multi-modal analysis by integrating transcriptomic and proteomic data from the same single cells. This integration reveals biological insights that would be impossible to obtain from either data type alone.

Data Integration Methods

Multi-Modal Integration Approaches



WNN Algorithm

For each cell, compute modality-specific weights:

Step 1: Calculate within-modality k-NN

- Find k nearest neighbors in RNA space
- Find k nearest neighbors in ADT space

Step 2: Compute modality weights

- Weight = consistency of neighborhood
- Higher weight if modality shows clear structure

Step 3: Construct weighted graph

- Combine RNA and ADT graphs
- Edge weight = $w_{\text{RNA}} \times d_{\text{RNA}} + w_{\text{ADT}} \times d_{\text{ADT}}$

Step 4: Clustering on WNN graph

CITE-seq Applications Across Biology

1 Immune Cell Profiling

- Resolve fine T cell subsets (CD4+, CD8+, Tregs, exhausted T cells)
- Characterize activation states (CD25, CD69, HLA-DR)
- Identify rare populations (double-negative T cells, $\gamma\delta$ T cells)
- Study immune responses in vaccination, infection, and cancer

Example: COVID-19 immune profiling, CAR-T cell monitoring

2 Tumor Microenvironment

- Map tumor-infiltrating lymphocytes (TILs)
- Identify immunosuppressive cells (MDSCs, Tregs, TAMs)
- Measure checkpoint molecules (PD-1, PD-L1, CTLA-4, TIM-3)
- Correlate immune infiltration with tumor genetics

Example: Predicting immunotherapy response, resistance mechanisms

3 Developmental Biology

- Track cell lineage markers during differentiation
- Validate developmental trajectories with surface markers
- Study stem cell heterogeneity (HSCs, MSCs, iPSCs)
- Measure differentiation stage markers (CD34, CD38, CD117)

Example: Hematopoiesis mapping, organoid development

4 Autoimmune & Inflammatory Disease

- Profile immune dysregulation in disease vs. healthy tissue
- Identify pathogenic cell states
- Track inflammatory markers (TNF- α , IL-6, IFN- γ receptors)
- Monitor treatment response over time

Example: Rheumatoid arthritis, inflammatory bowel disease, lupus

5 Neuroscience

- Classify neuronal subtypes using surface receptors
- Study neurotransmitter receptor expression
- Profile microglia activation states
- Investigate neuroinflammation in disease

Example: Alzheimer's disease, neurodegeneration, brain development

6 Sample Multiplexing (Cell Hashing)

- Pool multiple samples in one 10X run
- Use hashtag oligos (HTOs) for sample identification
- Reduce batch effects and costs
- Detect and remove doublets computationally

Example: Patient cohorts, treatment time courses, genetic perturbations

Advantages Over RNA-seq Alone

- **Resolving Ambiguous Clusters:** Protein markers can distinguish cell types that appear similar at the transcriptome level. For example, CD4+ and CD8+ T cells may cluster together based on RNA but are cleanly separated by CD4/CD8 protein expression.


- **Validation of Annotations:** Surface protein measurements provide independent validation of RNA-based cell type assignments, increasing confidence in biological interpretations.
- **Functional Markers:** Many activation markers (PD-1, Ki-67, phospho-proteins) are post-translationally regulated and not reliably detected at the RNA level.
- **Low-Dropout Measurements:** Protein counts are less sparse than RNA counts, providing more robust quantification of key markers that may be subject to dropout in scRNA-seq.

Computational Tools for CITE-seq Analysis

- **Seurat (R):** Integrated workflow with WNN analysis, multi-modal clustering, and visualization. Most widely used for CITE-seq data.
- **MUON (Python):** Multi-modal omics analysis framework with support for CITE-seq integration.
- **CiteFuse (R):** Specialized package for CITE-seq analysis including doublet detection and data normalization.
- **totalVI (Python/scvi-tools):** Probabilistic model for joint analysis of RNA and protein data with batch correction.
- **Azimuth:** Reference-based cell type annotation using CITE-seq reference atlases.

Future Directions and Extensions

CITE-seq technology continues to evolve with new applications and extensions including REAP-seq (RNA expression and protein sequencing), ASAP-seq (ATAC-seq with protein), TEA-seq (Transcriptome, Epitope, and Accessibility sequencing), and DOGMA-seq (DNA, Oligo, Gene expression, Methylation, and Accessibility sequencing). These multi-omic approaches promise even deeper insights into cellular identity and function by layering additional modalities onto single-cell measurements.

 CITE-seq represents a paradigm shift in single-cell biology by enabling simultaneous measurement of transcriptome and proteome, providing unprecedented resolution for cell type identification, functional state characterization, and disease profiling.