

Targeted Gene Panels

Overview

- Sequence specific set of genes related to condition
- Highly focused - typically 10-500 genes
- Very high coverage for selected regions (>500X)

Common Panel Types

Cancer

50-500 genes

Oncology hotspots

Cardio

50-200 genes

Heart conditions

Neuro

100-300 genes

Epilepsy, ataxia

Advantages

- Cost-effective (\$100-300)
- Very high depth
- Faster turnaround
- Detect low-frequency variants

Use Cases

- Hereditary cancer screening
- Pharmacogenetic testing
- Carrier screening
- Targeted diagnostics

Best for known genes associated with specific phenotypes

Principle: Selective Enrichment



- Target enrichment focuses sequencing on specific genomic regions of interest
- Reduces sequencing cost by 10-1000x compared to whole genome sequencing
- Increases depth of coverage for better variant detection
- Enables detection of low-frequency somatic variants (as low as 1-5%)

Capture Methods

1. Hybridization Capture (Solution-based)

Custom oligonucleotide probes (baits) complementary to target regions are mixed with fragmented DNA library. Target fragments hybridize to biotinylated probes and are captured using streptavidin-coated magnetic beads. Non-target DNA is washed away.

Examples: Agilent SureSelect, IDT xGen, Twist Bioscience

Best for: Larger panels (>100 genes), exome sequencing

2. Amplicon-based Sequencing (PCR)

Multiple primer pairs designed to amplify specific target regions simultaneously in a single multiplex PCR reaction. Amplified products are pooled and sequenced directly.

Examples: Illumina AmpliSeq, Ion Torrent AmpliSeq

Best for: Small-medium panels (10-200 genes), hotspot regions

3. Molecular Inversion Probes (MIPs)

Single-stranded DNA probes with sequences complementary to regions flanking the target. After hybridization, the probe circularizes around the target sequence, which is then amplified.

Best for: SNP genotyping, copy number variation detection

Targeted Panel Sequencing Workflow

1

DNA Extraction & QC

Extract high-quality genomic DNA from sample (blood, tissue, saliva). Assess quantity (10-500 ng typically required) and quality (DIN/RIN score).


2

Library Preparation

Fragment DNA to optimal size (150-300 bp). Attach adapters and unique molecular identifiers (UMIs) to enable sequencing and reduce PCR duplicates.

3

Target Enrichment



Method Comparison

Feature	Hybridization Capture	Amplicon-based
Input DNA	50-500 ng	10-50 ng
Uniformity	Excellent across targets	Variable (primer efficiency)
Target Size	Best for large panels (> 1 Mb)	Best for small panels (<500 kb)
Workflow Time	2-3 days	1 day

Feature	Hybridization Capture	Amplicon-based
Sensitivity	5-10% allele frequency	1-5% allele frequency
Cost per Sample	\$150-400	\$100-250
Best Application	Hereditary disease panels, exomes	Cancer hotspots, pharmacogenetics

Key Performance Metrics

Coverage Metrics

- **Mean coverage depth:** >500X
- **Target coverage:** >95% at 100X
- **Uniformity:** >80% bases within 0.2x mean
- **On-target rate:** >50% reads

Variant Detection

- **SNV sensitivity:** >99%
- **Indel sensitivity:** >95%
- **CNV detection:** Exon-level resolution
- **Somatic VAF:** As low as 1-5%

High depth sequencing enables confident detection of both germline and somatic variants with clinical-grade accuracy

Clinical Applications

- **Oncology:**Somatic mutation profiling for targeted therapy selection (e.g., EGFR, KRAS, BRAF in solid tumors)
- **Hereditary Cancer:**BRCA1/2, Lynch syndrome genes (MLH1, MSH2, MSH6, PMS2), Li-Fraumeni (TP53)
- **Cardiovascular:**Cardiomyopathy genes (MYH7, MYBPC3, TTN), arrhythmia panels (SCN5A, KCNQ1)
- **Neurology:**Epilepsy genes (SCN1A, KCNQ2), intellectual disability panels, muscular dystrophy genes

- **Pharmacogenomics:** Drug metabolism genes (CYP2D6, CYP2C19, TPMT, SLCO1B1) for personalized medication