

SAM/BAM Formats

SAM (Sequence Alignment/Map) Format

```
Header : @HD VN:1.6 S0:coordinate
@SQ SN:chr1 LN:248956422
Alignment: READ1 99 chr1 10001 60 76M = 10052 127 ACGT... I I I I...
```

SAM (Text)

- Human-readable
- Tab-delimited
- Large file size
- 11 mandatory fields

BAM (Binary)

- Compressed SAM
- ~3-5x smaller
- Faster to process
- Requires indexing (.bai)

Key SAM Fields

- QNAME - Read name
- FLAG - Bitwise flag (paired, mapped, reverse, etc.)
- RNAME - Reference sequence name (chromosome)
- POS - Alignment position
- MAPQ - Mapping quality score
- CIGAR - Alignment string (M=match, I=insertion, D=deletion)

1. QNAME (Query Name) - Read Name

QNAME: HWUSI-EAS100R:6:73:941:1973

QNAME: SRR123456.1

QNAME: READ_00001

Description: Unique identifier for each read generated by the sequencing instrument.

- Typically contains sequencer information and coordinate data
- Paired reads from the same DNA fragment share the same QNAME
- Format: InstrumentID:RunNumber:FlowCell:Lane:Tile:X-coord:Y-coord

2. FLAG - Bitwise Flag

FLAG Value Calculation Example

FLAG = 99 (0x63)
= 1 (paired) + 2 (properly paired) + 32 (mate reverse strand) + 64
(first in pair)

0x1 (1) = paired

0x2 (2) = properly paired

0x4 (4) = unmapped

0x8 (8) = mate unmapped

0x10 (16) = reverse strand

0x20 (32) = mate reverse

0x40 (64) = first in pair

0x80 (128) = second in pair

Description: Integer value representing various read attributes in binary format.


- Efficiently stores multiple attributes in a single number
- Each bit has independent meaning (verified using AND operation)
- Example: FLAG=99 → Paired read, properly mapped, first read, mate is reverse strand

3. RNAME & POS - Reference Sequence Name and Position

RNAME: chr1, chr2, chrX, chrY, chrM

POS: 10001 (1-based coordinate)

Chromosome Position Visualization



← chr1 start (1) Read mapping position chr1 end (248,956,422) →

Description:

- **RNAME:** Chromosome or contig name in reference genome (e.g., chr1, chr2, chrX)
- **POS:** 1-based coordinate where the read is mapped (chromosome starts at 1)
- '*' or '0' indicates an unmapped read
- BAM file indexing (.bai) enables rapid retrieval of specific regions

4. MAPQ - Mapping Quality Score

MAPQ Value	Meaning	Error Probability	Confidence
60	Very high quality	0.0001% (1/1,000,000)	✓✓✓✓✓
40	High quality	0.01% (1/10,000)	✓✓✓✓
20	Medium quality	1% (1/100)	✓✓✓
10	Low quality	10% (1/10)	✓✓
0	No mapping confidence	-	✗

Description: Phred-scaled mapping quality score ($-10 \times \log_{10}(P)$)

- Represents the probability that the read is mapped to an incorrect position
- MAPQ ≥ 30 : Generally considered reliable mapping
- MAPQ 0: Multi-mapping (maps equally well to multiple locations)
- Frequently used as filtering criteria in variant calling

💡 SAM/BAM formats are the standard formats for NGS data analysis, serving as the foundation for all sequencing analyses including variant analysis, RNA-seq, and ChIP-seq.

📌 5. CIGAR - Alignment String

CIGAR: 50M2I48M1D25M

50M

2I

48M

1D

25M

Visual Representation:

|||||

NM:i:0 MD:Z:76 AS:i:76

Field	Value	Description
QNAME	READ1	Read identifier
FLAG	99	Paired, properly mapped, first read, mate reverse strand
RNAME	chr1	Chromosome 1
POS	10001	Start position
MAPQ	60	Very high mapping quality
CIGAR	76M	76 bases perfectly aligned
RNEXT	=	Mate also mapped to same chromosome
PNEXT	10052	Mate start position
TLEN	127	Template length (insert size)
SEQ	ACGT...	Read sequence
QUAL	IIII...	Per-base quality scores (Phred+33)

7. BAM File Operations Example

```
# Convert SAM to BAM
samtools view -bS input.sam > output.bam

# Sort BAM
```

```
samtools sort output.bam -o sorted.bam
```

```
# Index BAM (creates .bai file)
```

```
samtools index sorted.bam
```

```
# Extract specific region
```

```
samtools view sorted.bam chr1:10000-20000
```

```
# View statistics
```

```
samtools flagstat sorted.bam
```

Description: BAM files are essential for efficient data processing.

- Only sorted BAM files can be indexed
- .bai index file enables rapid retrieval of specific regions
- Visualization tools like IGV require BAM + index
- Typical compression ratio is 1/3 to 1/5 compared to SAM