

# Library Preparation



## 1. DNA Fragmentation

Break DNA into smaller fragments (200-600 bp)



## 2. End Repair

Create blunt ends on DNA fragments



## 3. Adapter Ligation

Attach sequencing adapters to both ends



## 4. Size Selection

Select fragments of desired length



## 5. PCR Amplification

Amplify library for sequencing

### Quality Control

Check library size distribution and concentration using Bioanalyzer or TapeStation

### Critical Factors

Input DNA quality, fragmentation method, and adapter ligation efficiency