

Lipidomics: Comprehensive Analysis Guide

Lipid Classes

- Glycerophospholipids, sphingolipids
- Triacylglycerols, cholesterol esters
- 1000s of lipid species

Extraction Protocols

- Bligh-Dyer, Folch methods
- Biphasic extraction
- Lipid class-specific protocols

Separation Strategies

- Direct infusion (shotgun lipidomics)
- LC-MS with C8/C18 columns
- Supercritical fluid chromatography

Nomenclature

- Lipid MAPS classification
- Fatty acid composition notation
- Standardized reporting

1. Lipid Classes: Diversity and Structure

Lipids are a diverse group of hydrophobic or amphipathic molecules that play crucial roles in cellular structure, energy storage, and signaling. The mammalian lipidome comprises over 1,000 distinct lipid species across multiple classes.

Major Lipid Categories

Glycerophospholipids

Phosphatidylcholine (PC)
Phosphatidylethanolamine (PE)
Phosphatidylserine (PS)
Phosphatidylinositol (PI)

Sphingolipids

Ceramides (Cer)
Sphingomyelins (SM)
Glycosphingolipids
Gangliosides


Neutral Lipids

Triacylglycerols (TAG)
Diacylglycerols (DAG)
Cholesterol esters (CE)
Free fatty acids (FFA)

Sterols

Cholesterol
Oxysterols
Steroid hormones
Bile acids

Structural Representation: Glycerophospholipid Architecture



Example: Common Lipid Species

Lipid Class	Example Species	Notation	Biological Role
-------------	-----------------	----------	-----------------

Phosphatidylcholine	PC 16:0/18:1	PC(34:1)	Membrane structure, signaling
Ceramide	Cer d18:1/16:0	Cer(34:1)	Apoptosis, cell differentiation
Triacylglycerol	TAG 16:0/18:1/18:2	TAG(52:3)	Energy storage
Cholesterol Ester	CE 18:2	CE(18:2)	Cholesterol storage, transport

2. Extraction Protocols: Optimizing Lipid Recovery

Lipid extraction is the critical first step in lipidomics analysis. The choice of extraction method depends on the sample type, lipid classes of interest, and downstream analytical platform.

Classical Extraction Methods

Folch Method (1957)

Chloroform:Methanol (2:1, v/v)
Best for: Total lipid extraction
Recovery: >95% for most lipids
Volume ratio: 20:1 (solvent:sample)

Bligh-Dyer Method (1959)

Chloroform:Methanol:Water (1:2:0.8)
Best for: Aqueous samples
Lower solvent volume
Modified for tissue/cells

MTBE Method



Methyl-tert-butyl ether based
Less toxic than chloroform
Upper phase contains lipids
Excellent for neutral lipids

Butanol-Methanol

Alternative non-chlorinated
Good phospholipid recovery

Single-phase extraction
LC-MS compatible

Biphasic Extraction Process Flow



- **Temperature Control:** Keep samples at 4°C during extraction to prevent oxidation
- **Internal Standards:** Add deuterated standards before extraction for quantification
- **Sample Preparation:** Homogenize tissues thoroughly for reproducible extraction
- **Antioxidants:** Add butylated hydroxytoluene (BHT) to prevent lipid peroxidation
- **Storage:** Store lipid extracts at -80°C under nitrogen or argon

3. Separation Strategies: Analytical Approaches

Modern lipidomics employs various separation techniques coupled with mass spectrometry to achieve comprehensive lipid profiling. Each approach offers unique advantages for different analytical goals.

Comparison of Major Separation Techniques

Method	Principle	Advantages	Limitations	Best Application
Shotgun Lipidomics	Direct infusion ESI-MS/MS	High throughput Minimal sample prep Quantitative	Ion suppression Limited isomer separation	High-throughput screening Targeted analysis
LC-MS	Reverse-phase chromatography	Excellent resolution Reduces ion suppression Isomer separation	Longer analysis time Method development	Comprehensive profiling Complex mixtures

SFC-MS


Supercritical CO₂

Fast separation
Lipid class
separation
High efficiency

Specialized
equipment
Method optimization

Neutral lipid analysis
Isomer separation

LC-MS Workflow for Lipidomics



Column Selection for LC-MS Lipidomics

C8 Columns

Particle size: 1.7-3 µm
Length: 50-150 mm

C18 Columns

Particle size: 1.7-2.1 µm
Length: 100-250 mm

HILIC Columns

Phase: Hydrophilic interaction
Best for: Polar lipids (gangliosides,

Best for: Complex lipid mixtures, shorter analysis time
Retention: Moderate hydrophobic interaction

Best for: Maximum resolution, isomer separation
Retention: Strong hydrophobic interaction

cardiolipins)
Mobile phase: High organic content
Retention: Based on polarity

Typical LC-MS Method Parameters

Column: C18 (2.1 × 100 mm, 1.7 μm)

Mobile Phase A: Water:Acetonitrile (60:40) + 10 mM ammonium formate

Mobile Phase B: Isopropanol:Acetonitrile (90:10) + 10 mM ammonium formate

Flow Rate: 0.4 mL/min

Column Temperature: 55°C

Gradient: 40-100% B over 10 minutes


Ionization: ESI positive and negative mode

Mass Range: m/z 200-2000

4. Nomenclature: Standardized Lipid Annotation

Standardized nomenclature is essential for reproducible lipidomics research. The LIPID MAPS consortium has established a comprehensive classification system that enables consistent lipid identification and reporting across laboratories.

LIPID MAPS Classification System



Fatty Acid Composition Notation Examples

PC 34:1

Level 1 (Sum Composition): Total carbons:total double bonds

Example: Could be 16:0/18:1 or 16:1/18:0

PC 16:0_18:1

Level 2 (Molecular Species): Individual fatty acid compositions

sn position unknown

PC 16:0/18:1(9Z)

Level 3 (sn Position & Double Bond): Complete structural information

16:0 at sn-1, 18:1 with double bond at position 9 (Z configuration) at sn-2

PC O-16:0/18:1

Special Notation (Ether Lipid): O- indicates ether linkage

P- indicates plasmalogen (vinyl ether)

Common Lipid Abbreviations

Abbreviation	Full Name	LIPID MAPS Category	Example Species
PC	Phosphatidylcholine	Glycerophospholipids (GP)	PC 16:0/18:1, PC O-18:0/20:4
PE	Phosphatidylethanolamine	Glycerophospholipids (GP)	PE 18:0/20:4, PE P-16:0/22:6
SM	Sphingomyelin	Sphingolipids (SP)	SM d18:1/16:0, SM d18:1/24:1
Cer	Ceramide	Sphingolipids (SP)	Cer d18:1/16:0, Cer d18:1/24:0
TAG	Triacylglycerol	Glycerolipids (GL)	TAG 16:0/18:1/18:2, TAG 52:3
DAG	Diacylglycerol	Glycerolipids (GL)	DAG 18:1/18:1, DAG 36:2
CE	Cholesteryl Ester	Sterol Lipids (ST)	CE 18:2, CE 20:4
LPC	Lysophosphatidylcholine	Glycerophospholipids (GP)	LPC 18:0, LPC 20:4

Best Practices for Lipid Reporting

- **Use Standard Abbreviations:** Follow LIPID MAPS nomenclature for consistency
- **Specify Identification Level:** Indicate whether structure is fully resolved or putative

- **Report Mass Accuracy:** Include m/z values and mass errors (typically <5 ppm)
- **Document Retention Time:** Aids in identification verification and cross-study comparison
- **Include Internal Standards:** Report which standards were used for quantification
- **Provide Method Details:** Specify ionization mode, fragmentation patterns used for ID
- **Use LIPID MAPS ID:** Cross-reference with LIPID MAPS database when possible
- **Specify Isomer Information:** Note if cis/trans or sn-position was determined

Resources for Lipid Identification

LIPID MAPS

Comprehensive database
>47,000 lipid structures
Tools for structure drawing
MS/MS spectral libraries

LipidBlast

In silico MS/MS library
>200,000 spectra
Multiple ionization modes
Free download available

SwissLipids

Curated database
>777,000 lipid species
Extensive isomer coverage
Prediction tools

LipidSearch

Commercial software
Automated identification
Quantification tools
Statistical analysis

Key Takeaways for Successful Lipidomics Analysis

Comprehensive lipidomics requires careful attention to all analytical stages: from selecting appropriate extraction protocols for your sample type, to choosing optimal separation strategies for your lipid classes of interest, and finally to applying standardized nomenclature for reproducible reporting.

Integration of these methodological components ensures robust and meaningful biological insights.